



GENESIS

Genesis lab for diagnostic reagents

PRODUCT CATALOGUE 2026





CONTACT INFORMATION



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632 Fax : (+202) 44891632




info@genesis-egy.com



www.genesis-egy.com

Our Website scan Me:-





GENESIS CERTIFICATES



CMC Medical device, C/Horacio
Lengo n18 C.P29006, Malaga (Spain)



GENESIS



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com

PAGE 2



COMPANY PROFILE

GENESIS is an Egyptian Company Started up in 2014 for Biotechnology develops and manufactures high quality in diagnostic reagents.

OUR VALUE

is absolutely self-directed, with the influence of our Diagnostics expert team, violently researching team.

Giving – up our liquid chemistry reagents the maximum affordable stability which enhances linearity, reproducibility, reliability and sensitivity.

OUR VISION

To gain the respect and trust of our customers worldwide by continuously finding better solutions across all features of our products as one of the leading manufacturers of Diagnostics Products in the world.

OUR MISSION

provides the value of uniqueness.

We improve existing and add new qualities to products at every stage of their production, Delivers exclusive and supreme services stating partnership prospective ,supporting customers satisfaction guarantees successful over the world, commitment to quality.

PHILOSOPHY

Our philosophy in words:



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com

PAGE 3



GENESIS PRODUCTS

- 1. Clinical Chemistry

05

1.1. Substrates

1.2. Enzymes

1.3. Electrolyte

- 2. Hematology

13

- 3. Coagulation

14

- 4. Serology (Latex)

15

- 5. Blood Group

16

- 6. Bacteriology

18

- 7. TurbiLatex

19



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com




www.genesis-egy.com

PAGE 4

1. CLINICAL CHEMISTRY

1.1. SUBSTRATES



Albumin (BCG METHOD)

- Size : 2*100 & 4*100 & 2*50
- Method: End Point
- Wavelength: 578 nm
- Linearity: 7.0 g/dL
- Shelf life: 21 months
- Storage : 2-8°C



Bilirubin JENDRASSIK (Total, Direct)

- Size : 265 ml
- Method: End Point
- Wavelength:
Total : 578 nm & Direct : 546 nm
- Linearity: 20 mg/dL
- Shelf life: 24 months
- Storage : Store at 20-25°C



Bilirubin DMSO (Total, Direct)

- Size : 2*100 & 2*60 & 2*25
- Method: End Point
- Wavelength: 546 nm
- Linearity: 18 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C



Cholesterol (CHOD/PAP Method)

- Size : 4*25 & 2*50 & 4*50 & 2*25
- Method: End Point
- Wavelength: 546 nm
- Linearity: 1000 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C




Creatinine (Jaffe Method)

- Size : 2*100 & 4*100 & 2*50
- Method: Fixed-time
- Wavelength: 492 nm
- Linearity: 20 mg/dL
- Shelf life: 21 months
- Storage : 20-25°C




Glucose (GOD/PAP Method)

- Size : 2*100 & 4*100 & 2*50 & 4*250
- Method: End Point
- Wavelength: 546 nm
- Linearity: 500 mg/dL
- Shelf life: 24 months
- Storage : 2-8°C



HDL Cholesterol (Precipitating Method)


- Size : 1*50
- Method: End Point
- Wavelength: 546 nm
- Linearity: 150 mg/dL
- Shelf life: 24 months
- Storage : 15-25°C



Total Protein (BIURET Method)


- Size : 2*100 & 4*100 & 2*50
- Method: End Point
- Wavelength: 546 nm
- Linearity: 15.0 g/dL
- Shelf life: 21 months
- Storage : 2-8°C





Triglycerides (GPO / PAP Method)

- Size : 4*25 & 2*50 & 4*50 & 2*25
- Method: End-Point
- Wavelength: 546 nm
- Linearity: 600 mg/dL
- Shelf life: 21 months
- Storage : 20-25°C



Urea (Colorimetric Method)

- Size : 2*100 & 2*50
- Method: End Point
- Wavelength: 578 nm
- Linearity: 300 mg/dL
- Shelf life: 24 months
- Storage : 2-8°C



Urea UV (4+1)

- Size : 5*10 & 10*10 & 5*20 & 10*20 & 2*25
- Method: Fixed-time
- Wavelength: 340 nm
- Linearity: 300 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C



Uric Acid (Urease/PAP Method)

- Size : 4*25 & 2*50 & 4*50 & 2*25
- Method: End Point
- Wavelength: 546 nm
- Linearity: 25 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C



LDL (Direct Method)

- Size : 1*40 & 2*40
- Method: Fixed-time
- Wavelength: 578 nm
- Linearity: 400 mg/dl
- Shelf life: 20 months
- Storage : 2-8°C



HDL (Direct Method)

- Size : 1*40 & 2*40
- Method: Fixed-time
- Wavelength: 578 nm
- Linearity: 200 mg/dL
- Shelf life: 20 months
- Storage : 2-8°C



Micro total protein (M-TP)

- Size : 2*50 & 2*100 & 4*100
- Method: End point
- Wavelength: 578 nm
- Linearity: 400 mg/dL
- Shelf life: 24 months
- Storage : 2-8°C



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com

1.2. ENZYMES



Alkaline-Phosphatase (4+1) (DGKC Method)

- Size : 5*10 & 5*20 & 10*10 & 2*25
- Method: KINETIC
- Wavelength: 405 nm
- Linearity: 900 U/L
- Shelf life: 18 months
- Storage : 2-8°C



CK-NAC (4+1) (IFCC Method)

- Size : 5*5 & 10*5 & 1*25
- Method: KINETIC
- Wavelength: 340 nm
- Linearity: 400 U/L
- Shelf life: 21 months
- Storage : 2-8°C



CK-MB (4+1) (Immunoinhibition Method)

- Size : 5*5 & 10*5 & 1*25
- Method: Fixed-time
- Wavelength: 340 nm
- Linearity: 1000 U/L
- Shelf life: 18 months
- Storage : 2-8°C



GOT (AST) (4+1)

- Size : 5*10 & 10*10 & 5*20 & 10*20 & 2*25
- Method: KINETIC
- Wavelength: 340 nm
- Linearity: 400 U/L
- Shelf life: 21 months
- Storage : 2-8°C



GENESIS LAB FOR DIAGNOSTIC REAGENTS

1st industrial area, Obour City, Cairo, Egypt

(+202) 44891632 Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com



GPT (ALT) (4+1)

- Size : 5*10 & 10*10 & 5*20 & 10*20 & 2*25
- Method: KINETIC
- Wavelength: 340 nm
- Linearity: 400 U/L
- Shelf life: 21 months
- Storage : 2-8°C



LDH (4+1) (Pyruvate-Lactate Method)

- Size : 2*25 & 5*20
- Method: KINETIC
- Wavelength: 340 nm
- Linearity: 1250 U/L
- Shelf life: 21 months
- Storage : 2-8°C



α- Amylase (CNPG3)


- Size : 2*25 & 4*25
- Method: KINETIC
- Wavelength: 405 nm
- Linearity: 1500 U/L
- Shelf life: 18 months
- Storage : 2-8°C



γ- GT (IFCC Method)

- Size : 2*25 & 5*20
- Method: KINETIC
- Wavelength: 405 nm
- Linearity: 300 U/L
- Shelf life: 21 months
- Storage : 2-8°C

1.3. ELECTROLYTE



Calcium (Two Reagent) (OCPC Method)

- Size : 2*25 & 4*25 & 2*50
- Method: End-Point
- Wavelength: 546 nm
- Linearity: 20 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C




Calcium (Single Reagent)

- Size : 2*25 & 4*25 & 2*50
- Method: End-Point
- Wavelength: 620 nm
- Linearity: 20 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C



Chloride (Single Reagent)

- Size : 2*25 & 4*25
- Method: End point
- Wavelength: 500 nm
- Linearity: 190 mmol/L
- Shelf life: 21 months
- Storage : 2-8°C



Magnesium (Xylidyl Blue)

- Size : 2*25 & 4*25
- Method: End point
- Wavelength: 546 nm
- Linearity: 6 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C





Phosphorus UV

- Size : 2*25 & 4*25 & 2*50
- Method: End-Point
- Wavelength: 340 nm
- Linearity: 20 mg/dL
- Shelf life: 21 months
- Storage : 2-8°C



Potassium (Single Reagent)

- Size : 2*25 & 4*25
- Method: End-Point
- Wavelength: 630 nm
- Linearity: 20mmol/L
- Shelf life: 21 months
- Storage : 2-8°C



Sodium (Single Reagent)

- Size : 2*25 & 4*25
- Method: End point
- Wavelength: 630 nm
- Linearity: 180 mmol/L
- Shelf life: 21 months
- Storage : 2-8°C



Zinc

- Size : 2*25 & 4*25
- Method: End point
- Wavelength: 546 nm
- Linearity: 500 mg/dL
- Shelf life: 24 months
- Storage : 2-8°C



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632




info@genesis-egy.com




www.genesis-egy.com

2. HEMATOLOGY



**Hemoglobin
(Drabkin Reagent (9+1))**

- Size : 2*20 & 4*20
- Method: End-Point
- Wavelength: 546 nm
- Linearity: 20.0 g/dl
- Shelf life: 24 months
- Storage : 20-25°C



**Hemoglobin
(Drabkin Reagent (49+1))**

- Size : 2*20 & 4*20
- Method: End-Point
- Wavelength: 546 nm
- Linearity: 20.0 g/dl
- Shelf life: 24 months
- Storage : 20-25°C



**Hemoglobin
(Ready For Use)**

- Size : 2*250 & 4*250
- Method: End-Point
- Wavelength: 546 nm
- Linearity: 20.0 g/dl
- Shelf life: 24 months
- Storage : 20-25°C



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com





www.genesis-egy.com

3. COAGULATION



Prothrombin Time Reagent (PT)

- Size : 6*5 & 6*8 & 6*4
- Shelf life: 24 months
- Storage : 2-8 °C




Partial Thromboplastin Time Reagent (Aptt) (With Calcium Chloride)

- Size : 6*3 & 6*5
- Shelf life: 24 months
- Storage : 2-8°C




4. SEROLOGY (LATEX)



ASOT Latex

- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
- Shelf life: 30 months
- Storage : 2-8 °C




RF Latex

- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
- Shelf life: 30 months
- Storage : 2-8 °C




CRP Latex

- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
- Shelf life: 30 months
- Storage : 2-8 °C




5. BLOOD GROUP




Anti-A

- Size : 1*10 & 10*10
- Shelf life: 24 months
- Storage : 2-8 °C



Anti-B

- Size : 1*10 & 10*10
- Shelf life: 24 months
- Storage : 2-8 °C



Anti-D


- Size : 1*10 & 10*10
- Shelf life: 24 months
- Storage : 2-8 °C



Anti A,B,D

- Size : 3*10
- Shelf life: 24 months
- Storage : 2-8 °C





Anti-AB

- Size : 1*10
- Shelf life: 24 months
- Storage : 2-8 °C




Anti Human Globulin AHG


- Size : 1*10
- Shelf life: 24 months
- Storage : 2-8 °C

Bovine Serum Albumin (BSA) 22%

- Size : 1*10
- Shelf life: 24 months
- Storage : 2-8 °C



6. BACTERIOLOGY



Brucella Rose Bengal




- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
- Shelf life: 24 months
- Storage : 2-8 °C



Brucella A&M



- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
- Shelf life: 24 months
- Storage : 2-8 °C



Widal Test



- Size :
50 Test & 100 Test (With Controls)
50 Test & 100 Test (Without Controls)
(8 Antigen)
- Shelf life: 24 months
- Storage : 2-8 °C



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com

7. TURBILATEX



ASOT TurbiLatex



- Size : 50 Test & 100 Test & 200 Test
- Shelf life: 24 months
- Storage : 2-8 °C




CRP TurbiLatex




- Size : 50 Test & 100 Test & 200 Test
- Shelf life: 24 months
- Storage : 2-8 °C



RF TurbiLatex



- Size : 50 Test & 100 Test & 200 Test
- Shelf life: 24 months
- Storage : 2-8 °C



GENESIS LAB FOR DIAGNOSTIC REAGENTS



1st industrial area, Obour City, Cairo, Egypt



(+202) 44891632

Fax : (+202) 44891632



info@genesis-egy.com



www.genesis-egy.com



Ferritin TurbiLatex (9+1)

- Size : 25 Test & 50 Test & 100 Test & 200 Test
- Shelf life: 24 months
- Storage : 2-8 °C



Ferritin TurbiLatex (9+1) Table

Linear	Ferritin Conc	Linear	Ferritin Conc	Linear	Ferritin Conc
13.80	6	202.23	115	368.7	290
23.01	10	207.43	120	376.8	300
32.21	14	212.66	125	384.8	310
41.41	18	217.86	130	392.6	320
50.61	22	223.1	135	400.2	330
59.82	26	228.3	140	407.7	340
69.02	30	233.5	145	415.0	350
78.22	34	238.6	150	422.1	360
87.42	38	243.8	155	429.1	370
96.63	42	249.0	160	436.0	380
105.83	46	254.2	165	442.7	390
115.03	50	259.4	170	449.2	400
124.24	54	264.6	175	455.6	410
133.44	58	269.8	180	461.9	420
142.64	62	275.0	185	468.1	430
149.54	65	280.3	190	474.2	440
	68	285.5	195	480.1	450
	71	290.7	200	485.8	460
	74	295.9	205	491.4	470
	77	301.1	210	496.9	480
	80	306.3	215	502.3	490
	83	311.5	220	507.6	500
	86	316.7	225	512.8	510
	89	321.9	230	517.9	520



HbA1c

- Size : 50 Test & 100 Test & 200 Test
- Shelf life: 24 months
- Storage : 2-8 °C



HbA1c Table

X	A.C	X	A.C	X	A.C	X	A.C
52.9	3	152.3	6.6	255.4	10.2	261.1	10.8
54.6	3.1	155.2	6.7	257.8	10.3	263.3	10.9
56.4	3.2	158.1	6.8	259.8	10.4	265.8	11
58.1	3.3	161.0	6.9	262.0	10.5	268.3	11.1
59.8	3.4	163.9	7	264.2	10.6	270.8	11.2
61.5	3.5	166.8	7.1	266.4	10.7	273.3	11.3
63.2	3.6	169.7	7.2	268.6	10.8	275.8	11.4
64.9	3.7	172.6	7.3	270.8	10.9	278.3	11.5
66.6	3.8	175.5	7.4	273.0	11	280.8	11.6
68.3	3.9	178.4	7.5	275.2	11.1	283.3	11.7
70.0	4	181.3	7.6	277.4	11.2	285.8	11.8
71.7	4.1	184.2	7.7	279.6	11.3	288.3	11.9
73.4	4.2	187.1	7.8	281.8	11.4	290.8	12
75.1	4.3	190.0	7.9	284.0	11.5	293.3	12.1
76.8	4.4	192.9	8	286.2	11.6	295.8	12.2
78.5	4.5	195.8	8.1	288.4	11.7	298.3	12.3
80.2	4.6	198.7	8.2	290.6	11.8	300.8	12.4
81.9	4.7	201.6	8.3	292.8	11.9	303.3	12.5
83.6	4.8	204.5	8.4	295.0	12	305.8	12.6
85.3	4.9	207.4	8.5	297.2	12.1	308.3	12.7
87.0	5	210.3	8.6	299.4	12.2	310.8	12.8
88.7	5.1	213.2	8.7	301.6	12.3	313.3	12.9
90.4	5.2	216.1	8.8	303.8	12.4	315.8	13
92.1	5.3	219.0	8.9	306.0	12.5	318.3	13.1
93.8	5.4	221.9	9	308.2	12.6	320.8	13.2
95.5	5.5	224.8	9.1	310.4	12.7	323.3	13.3
97.2	5.6	227.7	9.2	312.6	12.8	325.8	13.4
98.9	5.7	230.6	9.3	314.8	12.9	328.3	13.5
100.6	5.8	233.5	9.4	317.0	13	330.8	13.6
102.3	5.9	236.4	9.5	319.2	13.1	333.3	13.7
104.0	6	239.3	9.6	321.4	13.2	335.8	13.8
105.7	6.1	242.2	9.7	323.6	13.3	338.3	13.9
107.4	6.2	245.1	9.8	325.8	13.4	340.8	14
109.1	6.3	248.0	9.9	328.0	13.5	343.3	14.1
110.8	6.4	250.9	10	330.2	13.6	345.8	14.2
112.5	6.5	253.8	10.1	332.4	13.7	348.3	14.3
114.2	6.6	256.7	10.2	334.6	13.8	350.8	14.4
115.9	6.7	259.6	10.3	336.8	13.9	353.3	14.5
117.6	6.8	262.5	10.4	339.0	14	355.8	14.6
119.3	6.9	265.4	10.5	341.2	14.1	358.3	14.7
121.0	7	268.3	10.6	343.4	14.2	360.8	14.8
122.7	7.1	271.2	10.7	345.6	14.3	363.3	14.9
124.4	7.2	274.1	10.8	347.8	14.4	365.8	15
126.1	7.3	277.0	10.9	350.0	14.5	368.3	15.1
127.8	7.4	280.0	11	352.2	14.6	370.8	15.2
129.5	7.5	282.9	11.1	354.4	14.7	373.3	15.3
131.2	7.6	285.8	11.2	356.6	14.8	375.8	15.4
132.9	7.7	288.7	11.3	358.8	14.9	378.3	15.5
134.6	7.8	291.6	11.4	361.0	15	380.8	15.6
136.3	7.9	294.5	11.5	363.2	15.1	383.3	15.7
138.0	8	297.4	11.6	365.4	15.2	385.8	15.8
139.7	8.1	300.3	11.7	367.6	15.3	388.3	15.9
141.4	8.2	303.2	11.8	369.8	15.4	390.8	16
143.1	8.3	306.1	11.9	372.0	15.5	393.3	16.1
144.8	8.4	309.0	12	374.2	15.6	395.8	16.2
146.5	8.5	311.9	12.1	376.4	15.7	398.3	16.3
148.2	8.6	314.8	12.2	378.6	15.8	400.8	16.4
150.0	8.7	317.7	12.3	380.8	15.9	403.3	16.5
151.7	8.8	320.6	12.4	383.0	16	405.8	16.6
153.4	8.9	323.5	12.5	385.2	16.1	408.3	16.7
155.1	9	326.4	12.6	387.4	16.2	410.8	16.8
156.8	9.1	329.3	12.7	389.6	16.3	413.3	16.9
158.5	9.2	332.2	12.8	391.8	16.4	415.8	17
160.2	9.3	335.1	12.9	394.0	16.5	418.3	17.1
161.9	9.4	338.0	13	396.2	16.6	420.8	17.2
163.6	9.5	340.9	13.1	398.4	16.7	423.3	17.3
165.3	9.6	343.8	13.2	400.6	16.8	425.8	17.4
167.0	9.7	346.7	13.3	402.8	16.9	428.3	17.5
168.7	9.8	349.6	13.4	405.0	17	430.8	17.6
170.4	9.9	352.5	13.5	407.2	17.1	433.3	17.7
172.1	10	355.4	13.6	409.4	17.2	435.8	17.8
173.8	10.1	358.3	13.7	411.6	17.3	438.3	17.9
175.5	10.2	361.2	13.8	413.8	17.4	440.8	18
177.2	10.3	364.1	13.9	416.0	17.5	443.3	18.1
178.9	10.4	367.0	14	418.2	17.6	445.8	18.2
180.6	10.5	369.9	14.1	420.4	17.7	448.3	18.3
182.3	10.6	372.8	14.2	422.6	17.8	450.8	18.4
184.0	10.7	375.7	14.3	424.8	17.9	453.3	18.5
185.7	10.8	378.6	14.4	427.0	18	455.8	18.6
187.4	10.9	381.5	14.5	429.2	18.1	458.3	18.7
189.1	11	384.4	14.6	431.4	18.2	460.8	18.8
190.8	11.1	387.3	14.7	433.6	18.3	463.3	18.9
192.5	11.2	390.2	14.8	435.8	18.4	465.8	19
194.2	11.3	393.1	14.9	438.0	18.5	468.3	19.1
195.9	11.4	396.0	15	440.2	18.6	470.8	19.2
197.6	11.5	398.9	15.1	442.4	18.7	473.3	19.3
199.3	11.6	401.8	15.2	444.6	18.8	475.8	19.4
201.0	11.7	404.7	15.3	446.8	18.9	478.3	19.5
202.7	11.8	407.6	15.4	449.0	19	480.8	19.6
204.4	11.9	410.5	15.5	451.2	19.1	483.3	19.7
206.1	12	413.4	15.6	453.4	19.2	485.8	19.8
207.8	12.1	416.3	15.7	455.6	19.3	488.3	19.9
209.5	12.2	419.2	15.8	457.8	19.4	490.8	20
211.2	12.3	422.1	15.9	460.0	19.5	493.3	20.1
212.9	12.4	425.0	16	462.2	19.6	495.8	20.2
214.6	12.5	427.9	16.1	464.4	19.7	498.3	20.3
216.3	12.6	430.8	16.2	466.6	19.8	500.8	20.4
218.0	12.7	433.7	16.3	468.8	19.9	503.3	20.5
219.7	12.8	436.6	16.4	471.0	20	505.8	20.6
221.4	12.9	439.5	16.5	473.2	20.1	508.3	20.7
223.1	13	442.4	16.6	475.4	20.2	510.8	20.8
224.8	13.1	445.3	16.7	477.6	20.3	513.3	20.9

Arab Republic of Egypt
Egyptian Drug Authority
C.A. of Operations

0199298

0199298

CERTIFICATE OF GMP COMPLIANCE

Certificate Number: M-1212/2025
Part 1:

EDA confirms the following:
Name of Manufacturing Site: Genesis Lab for Diagnostics
Address of Manufacturing Site: Rented Factory - Unit no. 356 in the 1st Floor - Talaat Harb Building on Parts no. (6,7,8,9)-block no. (12017)- Obour Industrial Complex - Obour city.
License Number: LVPD2022005/03/2022
Purpose of the certificate: To be introduced to Central Administration of Medical Devices
From the knowledge gained during inspection of this manufacturer, the latest of which was conducted in 16/01/2025 it is considered that it complies with the principles and guidelines of Good Manufacturing Practice according to WHO guidelines/ other approved international guidelines for the production lines listed in the part 2.

Part 2:

Production line

Diagnostic Kits:

1. Solutions (clinical chemistry)

- The responsibility for the quality of the individual batches the in vitro diagnostic (IVDs) manufactured through this process lies with the manufacturer.
- This certificate reflects the status of the manufacturing site at the time of the inspection noted above and the period of validity may be reduced or extended using regularity risk management principle by an entry in the Restrictions or Clarifying remarks field.
- The certificate remains valid until 15/01/2026 it becomes invalid if the activities certified here are changed or if the site is no longer considered to be in compliance with GMP.

Dr. Sally Abd El Rasoul
Manager of Administration of Medical Devices & In vitro diagnostic inspection

Dr. Omnia Abdallah
General manager of general Administration of factory inspection

Authenticated,
Dr. Amira Mahgoub
Head of Central Administration of Operations

Issue Date of certificate: 17/02/2025
The authenticity of this certificate can be verified via scanning the QR-Code at the upper right corner

CAO.037.01 Issue/Rev.No.3/1 Issue date: 12/11/2023 Rev. date: 03/09/2024

ABC CERTIFICATE

**ABC CERTIFICATION hereby certifies that
GENESIS LAB FOR DIAGNOSTIC REAGENTS**

UNIT NO. (356) / FIRST FLOOR OF TALAAT HARB BUILDING / OBOUR INDUSTRIAL COMPLEX PROJECT/ OBOUR CITY/ QALYUBIA - EGYPT.

**Has been audited by ABC CERTIFICATION and has implemented
Medical Devices - Quality Management Systems**

ISO 13485:2016

Scope of Certification
IVD DESIGN, DEVELOPMENT AND MANUFACTURING.

Code: A.1.4
Original Date: 20.11.2023
Re-certification: 19.11.2026
Approval Date: 15.11.2025
Valid Period : 20.11.2025 - 19.11.2026
Certificate No: ABC-281-MD23

IAF EGAC Accredited MDQMS Certification CAB # 0123109

ABC MDQMS Certification CAB # 0123109

OG&MM CHIEF EXECUTIVE

UK Office: Kemp House, 152-160 City Road, London EC1V 2NX, United Kingdom
Egypt Office: 29 Al-Obour Building, Salah Salem St, Floor 1, Flat 14, Cairo, Egypt
Tel.: +2 02 24015033 - 01555 0555 26 - +44 207 0199019/6
Email: info@abccerti.com Website: www.abccerti.com

CMC⁺ EU REP⁺
Medical Devices Spain

Hereby confirms that CMC Medical Devices & Drugs S.L. located in C/Horacio Lengo, 18, 29006, Málaga, Spain and with SRN ES-AR-000000293 is the Authorised Representative (EC REP) of

GENESIS LAB FOR DIAGNOSTIC REAGENTS
1st, industrial area, obour city, Egypt

The certificate remains valid until the expiration date of the EC REP agreement or until manufacturing conditions, the quality system or relevant legislation are modified.

The validity is subject to positive results of periodic surveillance audits.

The product liability rests with the manufacturer in accordance with the applicable directive/regulation mentioned in Annex I of this certificate and the applicable standards.

After fulfilling of the relevant EU legislation requirements, the manufacturer shall affix relevant CE marking to the relevant models of the medical devices mentioned below.

The table in the Annex I lists, where applicable, the registration status with the Spanish MOH (AEMPS) and their registration number (RPS).

CE

CMC⁺

Issue date: 18th of November 2025
Manuel Mateos
Chief Executive Officer

Expiration date: 17th of October 2026
Certificate #: ECREP20251118.1
Ver: CERT-EC-2407
cmcmedicaldevices.com

10 YEARS



“Quality, accuracy, and trust — at the core of everything we deliver.”

“Your trusted partner in diagnostic excellence.”

**THANK
YOU**